Package 'DBSverify'

August 1, 2021

Title Analyze Next Generation Sequencing Data to Verify Somatic Double-BaseSubstituion (DBS) Calls Description DBS verify stands for "Doublet Base Substitution verify". The main function is Read_DBS_VCF_and_BAMs_to_verify_DBSs. This reads a VCF (variant call format) file containing somatic DBS mutations and examines the supporting reads in the corresponding tumor and normal BAM files to assess wither each DBS is likely real, i.e. arose at one time on a single chromosome. This is as opposed to being two adjacent single base mutations on homologous chromosomes, or adjacent mutations that occurred on the same chromosome but at detectably different times. Many somatic mutation callers do not call DBSs at all, and those that do sometimes make obvious errors, such as calling DBSs that consist of a germline SNP next to a somatic single base mutation. Because of the way the functionality of this package is used, the main function for most users will be Read_DBS_VCF_and_BAMs_to_verify_DBSs. License GPL-3 **Encoding** UTF-8 LazyData true **Roxygen** list(markdown = TRUE) RoxygenNote 7.1.1 Imports stringr, ICAMS (>= 2.3.17), data.table Suggests spelling, testthat, dplyr, tibble Language en-US **R** topics documented: CategorizeReads

	DBS_conclusion_1_row	2
	GetAllBAMSlicesSamtools	3
	PCAWG_read_table_and_evaluate_DBS	4
	ReadSamfile	5
	Read_DBS_VCF_and_BAMs_to_verify_DBSs	6
	Read_SBS_VCF_and_BAMs_to_verify_DBSs	8
	Slice2ReadSupport	9
Index		10

CategorizeReads

Categorizes the reads in one sam file.

Description

Categorizes the reads in one sam file.

Usage

```
CategorizeReads(sam, POS, REF, ALT)
```

Arguments

sam	An in-memory version of a sam file as data.frame.
POS	The position in the sam at which the DBS starts.
REF	The reference allele.
ALT	The putative alternative allele.

Value

a named character vector, the names are the read names in the sam file. Each element is one of "Read supports only 1st position" "Read supports only 2nd position", "WT read" "Mut read"

@keywords internal

```
DBS_conclusion_1_row
```

Examine 1 row of a DBS-only VCF that was already processed by Slice2ReadSupport to decide if the DBS is likely real.

Description

Examine 1 row of a DBS-only VCF that was already processed by Slice2ReadSupport to decide if the DBS is likely real.

Usage

```
DBS_conclusion_1_row(row, germlineCutOff = 0.2, max.half.support.T.reads = 1)
```

GetAllBAMSlicesSamtools 3

Arguments

row

One row of a DBS VCF (as a data.frame) already processed by Slice2ReadSupport (so that the fields vcf\$NreadSupport and vcf\$TreadSupport are pop-

ulated).

germlineCutOff

If this proportion of normal reads show one or the other variant (or both variants), consider this a germline variant or partial germline variant.

max.half.support.T.reads

Do not tolerate more than this number of reads in the tumor that support one but not both mutated positions.

Value

A character string indicating the conclusion about the putative DBS.

GetAllBAMSlicesSamtools

Get and save all slices of BAM as specified by positions in a VCF table.

Description

Get and save all slices of BAM as specified by positions in a VCF table.

Usage

```
GetAllBAMSlicesSamtools(
  vcf,
  bam.name,
  padding = 10,
  where.to.put.slices = tempfile(),
  verbose = 0
)
```

Arguments

vcf A VCF (Variant Call Format) "file" as a data.table or similar.

bam. name The file name (path) to the BAM file to slice.

padding How many base pairs on each side of the first base of the DBS to to keep in the

BAM slices.

where.to.put.slices

If NULL, create a temporary directory to store the slices. Otherwise, a character string that specifies a directory in which to store the BAM slices. This is

directory is created if necessary.

verbose If > 0 print a message when starting the number of slices generated every verbose

slices.

Value

A character string that specifies the directory the containing the BAM slices. The slices are stored as SAM files.

```
PCAWG_read_table_and_evaluate_DBS
```

Read a table that specifies how to process VCFs and mini-BAMS to evaluate DBS calls, for PCAWG Collaboratory data.

Description

Read a table that specifies how to process VCFs and mini-BAMS to evaluate DBS calls, for PCAWG Collaboratory data.

Usage

```
PCAWG_read_table_and_evaluate_DBS(
   in.table,
   in.vcf.dir,
   minibam.dir,
   out.vcf.dir = in.vcf.dir,
   bam.suffix = "_dbs_srt",
   verbose = 1
)
```

Arguments

in.table	the file path of the table to process; in production,/DBSverify/data-raw/collaboratory_bams_full_20 for testing/DBSverify/data-raw/short_collaboratory_bams.csv.
in.vcf.dir	The path to the directory containing the DBS VCF files.
minibam.dir	The path to the directory containing the mini BAMs.
out.vcf.dir	The path to the directory in which to put the "evaluated" DBS VCF files.
bam.suffix	String to add to end of BAM file name; depends on the conventions used by the script (run on the Collaboratory) that generated the miniBAMs.
verbose	If > 0 generate some progress messages.

Details

This is a specialized function for processing PCAWG data from the ICGC (International Cancer Genome) "Collaboratory" cloud computing system, once the miniBAMs have been created in the Collaboratory and downloaded. The <code>in.table</code> and associated BED files were used to specify the contents of the miniBAMs. The result consists of the "evaluated" DBS VCF files. The naming of the input and output VCF files and the miniBAMs is governed by the contents of <code>in.table</code>, with the VCF file names incorporating the "aliquot_id" and the miniBAM names based on the <code>icgc_donor_id</code> and the T_Specimen ID and N_Specimen ID.

ReadSamfile 5

ReadSamfile Read a "SAM file", discarding some reads that cannot be interpreted for our purposes.

Description

Read a "SAM file", discarding some reads that cannot be interpreted for our purposes.

Usage

ReadSamfile (filename)

Arguments

filename The name of the SAM file to read.

Details

SAM stands for "Sequence Alignment Map", a text file that represents aligned next-generation sequencing reads (https://en.wikipedia.org/wiki/SAM_(file_format)). Only keep reads that satisfy certain conditions:

- Mate pair maps to same chromosome
- Mapping quality >= 30
- "FLAG" < 256
- The CIGAR string is only \d+M (one or more digits followed M, and nothing else before or after). This means there are no insertions or deletions in the read versus the reference and there is no soft clipping. This function cannot keep track of read locations after insertions or deletions.

Value

A a list with the elements:

- good.reads, data.frame with with column names for the first 11 columns as specified in https://en.wikipedia.org/wiki/SAM_(file_format) with one row per read. The result contains only the reads that meet certain conditions.
- total.depth, the initial depth including "bad" reads
- reads.with.bad.FLAG. These are reads with FLAG >= 256, which marks reads that (i) are "not primary alignment" (ii) failed vendor QC (iii) are PCR or optical duplicates (iv) are supplementary alignments (e.g. split, split/inverted read). See https://broadinstitute.github.io/picard/explain-flags.html
- reads.with.bad.CIGAR, reads with CIGAR string that indicates an indel in the read or clipping.
- reads.with.bad.MAPQ, reads with MAPQ < 30.
- reads.with.bad.Mate_CHROM, reads with a mate on a different chromosome.

```
Read_DBS_VCF_and_BAMs_to_verify_DBSs
```

Determine whether sequencing reads in fact support (candidate) DBSs present in a VCF file.

Description

Determine whether sequencing reads in fact support (candidate) DBSs present in a VCF file.

Usage

```
Read_DBS_VCF_and_BAMs_to_verify_DBSs(
   input.vcf,
   Nbam.name,
   Tbam.name,
   N.slice.dir = tempfile(),
   T.slice.dir = tempfile(),
   unlink.slice.dir = TRUE,
   exclude.SBSs = TRUE,
   verbose = 0,
   outfile = NULL,
   filter.status = "PASS"
)
```

filter.status.

Arguments

	input.vcf	If a character string, then the path to a VCF file; otherwise A a single VCF "file" as a data.frame or similar object.
	Nbam.name	The name of the BAM file for the normal sample corresponding to ${\tt vcf.name}.$
	Tbam.name	The name of the BAM file for the tumor sample corresponding to ${\tt vcf.name}.$
	N.slice.dir	Directory for the slices of the normal BAM. Created if necessary.
	T.slice.dir	Directory for the slices of the tumor BAM. Created if necessary. Must be different than ${\tt N.slice.dir}$.
	unlink.slice	.dir
		If TRUE unlink N.slice.dir and T.slice.dir before return.
	exclude.SBSs	If TRUE silently filter out (exclude) SBSs in the input VCF. This makes sense if the the VCF is from a caller (like Mutect or the Hartwig Medical Foundation caller) that calls both SBSs and DBS.
	verbose	If > 0 print a message when starting the number of slices generated every <code>verbose</code> slices.
	outfile	If not NULL then write the "evaluated" VCF to outfile; otherwise write it to paste0 (input.vcf (vcf.name, "_evaluated.vcf"). Must be non-NULL if input.vcf is not a file path.
filter.status		
		If not NULL only keep rows where the FILTER column in the VCF is equal to

Details

Creates a new VCF file. This VCF file has no data rows if there were no DBSs to analyze. Otherwise, this VCF contains some additional columns. Any SBSs or indels in the input are silently ignored, and no attempt is made to merge adjacent SBSs.

- 1. NreadSupport With regard to the two positions of the DBS in the normal BAM, a string with 4 numbers separated by ":", with the numbers indicating respectively:
 - the number of reads that are reference sequence at both positions of the DBS,
 - the number of reads that that have the alternative allele only at the 1st position of the DBS
 - the number of reads that have the alternative allele only at the second position of the DBS,
 and
 - the number of reads that have the alternative alleles at both positions of the DBS.
- 2. TreadSupport Information analogous to that in NreadSupport, for the tumor BAM.
- 3. num_bad_mapped_reads The total number of tumor reads with MAPQ < 30 or with a mate on a different chromosome. If there are many badly mapped reads in the slice the slice may represent a segmental duplication.
- 4. num_bad_mapped_DBS_reads The number of tumor reads with the putative DBS but with MAPQ < 30 or a mate on a different chromosome. If many badly mapped reads support DBSs the DBS might results from mismapped reads in a segmental duplication.
- 5. DBSconclusion A string that describes whether the DBSs is believable ("True DBS"), or if the DBS is not believable, a string that describes why not.

The decision in DBSconclusion is based on multiple criteria. I also suggest relying on any available upstream filtering of SBSs that get merged into DBSs, as well as upstream filtering of DBSs. It is difficult to capture all the possible characteristics of likely miscalled DBSs, especially if they stem from mismapped reads that nevertheless have high MAPQ (mapping quality). The code that implements these criteria is in DBS_conclusion_1_row, which depends on classification of individual reads in ReadSamfile. The filtering of reads in the input SAM files is described in the documentation for ReadSamfile.

Once the reads to analyze are selected, additional criteria include

- There must be \geq 5 normal reads at the site of the putative tumor DBS.
- The normal reads at each separate position of the DBS must have < 20% variant calls.
- < 2 normal reads have the DBS.
- At least 2 tumor reads have the DBS.
- There are more well-mapped tumor reads than badly mapped tumor reads at the site of the DBS.
- There are more well-mapped tumor reads than badly mapped tumor reads that contain the DBS.
- If 1 normal read supports the DBS, then there must be a statistically greater proportion of tumor reads supporting the DBS (by Fisher's test).

Value

Invisibly, a list with the elements

- 1. The name of the DBS-only VCF file created.
- 2. The in-memory representation of the DBS VCF as a data.table.
- 3. The name of the directory with the normal SAM slices, if unlink.slice.dir is FALSE.
- 4. The name of the directory with the tumor SAM slices, if unlink.slice.dir is FALSE.

```
Read_SBS_VCF_and_BAMs_to_verify_DBSs
```

Aggressively merge SBSs into DBSs then determine whether sequencing reads in fact support DBSs.

Description

Aggressively merge SBSs into DBSs then determine whether sequencing reads in fact support DBSs.

Usage

```
Read_SBS_VCF_and_BAMs_to_verify_DBSs(
   input.vcf,
   Nbam.name,
   Tbam.name,
   variant.caller,
   N.slice.dir = tempfile(),
   T.slice.dir = tempfile(),
   unlink.slice.dir = TRUE,
   outfile = NULL,
   verbose = 0
)
```

Arguments

input.vcf	If a character string, then the path to a VCF file; otherwise A a single VCF "file" as a data.frame or similar object.			
Nbam.name	The name of the BAM file for the normal sample corresponding to vcf.name.			
Tbam.name	The name of the BAM file for the tumor sample corresponding to vcf.name.			
variant.caller				
	One of "strelka", "PCAWG", or "unknown". Merging adjacent SBS is done by ReadAndSplitVCFs in the ICAMS package.			
N.slice.dir	Directory for the slices of the normal BAM. Created if necessary.			
T.slice.dir	Directory for the slices of the tumor BAM. Created if necessary. Must be different than N.slice.dir.			
unlink.slice.dir				
	If TRUE unlink N.slice.dir and T.slice.dir before return.			
outfile	If not NULL then write the "evaluated" VCF to outfile; otherwise write it to paste0 (input.vcf (vcf.name, "_evaluated.vcf"). Must be non-NULL if input.vcf is not a file path.			
verbose	If > 0 print a message when starting the number of slices generated every <code>verbose</code> slices.			

Details

Note: argument input.vcf must be a file path. This function creates a new VCF file. See Read_DBS_VCF_and_BAMs_to_verify_DBSs for details.

Slice2ReadSupport 9

Value

 ${\bf Same~as~Read_DBS_VCF_and_BAMs_to_verify_DBSs.}$

Slice2ReadSupport Calculate the support for a putative DBS from the SAM slice containing the overlapping reads.

Description

Calculate the support for a putative DBS from the SAM slice containing the overlapping reads.

Usage

```
Slice2ReadSupport(slice.dir, CHROM, POS, REF, ALT)
```

Arguments

slice.dir The directory containing the SAM slices.

CHROM The chromosome identifier.

POS The first position of the DBS.

REF The reference variant.

ALT The alternate variant.

Index